

GLASGOW'S MEM, 1X*WITH L-GLUTAMINE & SODIUM BICARBONATE, WITHOUT TRYPTOSE PHOSPHATE BROTH, STERILE***Catalogue no.:** 322-005-****Shipping:** Ambient**Lot no.:** 322005054**Storage:** Store at 2 to 8°C. Protected from**Expiration date:** 08/2025

light

MSDS: Available upon request

TEST	SPECIFICATION	UNITS	RESULTS
Appearance	Clear red solution	NA	Conform
pH	7.2 ± 0.2	NA	7.29
Endotoxin	≤ 0.25	EU/ml	Conform
Mycoplasma	Tested Negative	NA	Negative
Sterility by USP	Pass	NA	Pass
Osmolality	290 ± 30	mOsm/Kg H ₂ O	299

ADDITIONAL INFORMATION:

STABILITY: This liquid media has a data supporting a shelf life of 18 months when properly stored. Avoid freezing of this media, as it may cause precipitation upon thawing.

PREPARATION: Add the appropriate amount of supplements (serum, antibiotics, growth factors) to sterile media in order to obtain the required working concentration, and adjust pH. Mix well and sterile filter (optional, reduces chances of contamination). Working concentration: selection of the proper working concentration is a balance between strong antibiotic effect and the concern for negative effects to the cells of importance. The working concentration should be varied above or below the suggested value to optimize for the cells.